

Product Datasheet

Haptoglobin Antibody HM707

Product Name	Haptoglobin Antibody
Catalog Number	HM707
Species/Host	Mouse
Tested Applications	LF
Immunogen	Purified Human Haptoglobin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Human Haptoglobin
Matched Pairs	HM706 (capture)
Preservatives	0.09 % Sodium Azide