

Product Datasheet

Cancer Antigen 72-4 (CA72-4) Antibody HM711

Product Name	Cancer Antigen 72-4 (CA72-4) Antibody
Catalog Number	HM711
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Purified human CA 72-4
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	Phosphate Buffered Saline, pH 7.0
Specificity	Human CA 72-4
Matched Pairs	HM710 (detection)
Preservatives	0.09% Sodium Azide