

Product Datasheet

Norovirus genogroup II Antibody HM713

Product Name	Norovirus genogroup II Antibody
Catalog Number	HM713
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Norovirus GII.4 VLP
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Norovirus GII
Matched Pairs	HM712 (capture)
Preservatives	0.09% Sodium Azide