

Product Datasheet

Anti-Mullerian Hormone (AMH) Antibody HM740

Product Name	Anti-Mullerian Hormone (AMH) Antibody
Catalog Number	HM740
Species/Host	Mouse
Tested Applications	LF, ELISA, CLIA
Immunogen	Recombinant AMH
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	>95%
Buffer	Phosphate Buffered Saline, pH 7.0
Specificity	AMH
Matched Pairs	HM1050 (detection)
Preservatives	0.09% Sodium Azide