

Product Datasheet

Hemoglobin (Hb) Antibody HM753

Product Name	Hemoglobin (Hb) Antibody
Catalog Number	HM753
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Human hemoglobin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Human hemoglobin
Matched Pairs	HM061 or HM782 (detection)
Preservatives	0.09% Sodium Azide