

## Product Datasheet

### Gastrin 17 (G-17) Antibody HM774

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Gastrin 17 (G-17) Antibody  |
| <b>Catalog Number</b>      | HM774   |
| <b>Species/Host</b>        | Mouse   |
| <b>Tested Applications</b> | LF, ELISA, CLIA, IFA  |
| <b>Immunogen</b>           | Recombinant Gastrin 17  |
| <b>Brand</b>               | EastCoast Bio   |
| <b>Form/Appearance</b>     | Protein A Purified Monoclonal Antibody                                    |
| <b>Storage</b>             | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| <b>Clonality</b>           | Monoclonal  |
| <b>Purity</b>              | >95% by SDS-PAGE  |
| <b>Buffer</b>              | 10 mM Phosphate Buffered Saline, pH 7.4                                   |
| <b>Specificity</b>         | Gastrin 17  |
| <b>Matched Pairs</b>       | HM775 or HM776 (detection)  |
| <b>Preservatives</b>       | 0.02% Proclin 300   |