

## Product Datasheet

### Gastrin 17 (G-17) Antibody HM774

<b>Product Name</b>	Gastrin 17 (G-17) Antibody
<b>Catalog Number</b>	HM774
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA, CLIA, IFA
<b>Immunogen</b>	Recombinant Gastrin 17
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Specificity</b>	Gastrin 17
<b>Matched Pairs</b>	HM775 or HM776 (detection)
<b>Preservatives</b>	0.02% Proclin 300