

Product Datasheet

Feline Immunodeficiency Virus (FIV-p24) Antibody HM780

| | |
|----------------------------|--|
| Product Name | Feline Immunodeficiency Virus (FIV-p24) Antibody |
| Catalog Number | HM780 |
| Species/Host | Mouse |
| Tested Applications | ELISA, WB, IHC, IFA |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | 2-8°C |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.2 |
| Cross Reactivity | No reactivity with FeLV, Feline Herpesvirus type 1, Feline Coronavirus, Feline Calicivirus |
| Specificity | FIV p24 (mature protein) |
| Matched Pairs | HM781 (detection) |
| Preservatives | 0.05% Sodium Azide |