

Product Datasheet

Progastrin-releasing Peptide (ProGRP) Antibody HM785

Product Name	Progastrin-releasing Peptide (ProGRP) Antibody
Catalog Number	HM785
Species/Host	Mouse
Tested Applications	ELISA, CLIA, CMIA
Immunogen	Recombinant ProGRP
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffer with 150 mM NaCl, pH 7.4
Matched Pairs	HM784 (capture)
Preservatives	0.1% Proclin 300