

## Product Datasheet

### Human IL-17F ELISA Kit 850.980

<b>Product Name</b>	Human IL-17F ELISA Kit
<b>Catalog Number</b>	850.980
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	15.6-500 pg/ml
<b>Sensitivity</b>	6.0 pg/ml
<b>Tested Applications</b>	ELISA
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reaction with IL-17A, IL-17B, IL-17D, IL-17E, IL-5, IL-8 and Perforin. Recognizes IL-17A/F.
<b>Specificity</b>	Recognizes both natural and recombinant human IL-17F
<b>Sample Types</b>	Serum, Plasma, Cell culture supernatant