

## Product Datasheet

### Parvovirus, Feline Antibody HM790

<b>Product Name</b>	Parvovirus, Feline Antibody
<b>Catalog Number</b>	HM790
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, IFA
<b>Immunogen</b>	Feline Parvovirus
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Cross Reactivity</b>	Canine Parvovirus
<b>Specificity</b>	Feline Parvovirus
<b>Matched Pairs</b>	HM790 (self-pairs) or HM789 (capture)
<b>Preservatives</b>	0.05% Sodium Azide