

Product Datasheet

Powassan Virus (POWV-NS1) Antibody HM794

Product Name	Powassan Virus (POWV-NS1) Antibody
Catalog Number	HM794
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant Powassan Virus
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM791, HM792 or HM793 (detection)
Preservatives	0.05% Sodium Azide