

## Product Datasheet

### Astrovirus Antibody HM803

<b>Product Name</b>	Astrovirus Antibody
<b>Catalog Number</b>	HM803
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA
<b>Immunogen</b>	Astrovirus
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Cross Reactivity</b>	No cross reactivity with Adenovirus, C. difficile, Cryptosporidium, Enterovirus, E. coli, Giardia, H. pylori, L. monocytogenes, Norovirus, Rotavirus, Y. enterocolitica
<b>Matched Pairs</b>	HM804 (capture or detection)
<b>Preservatives</b>	0.09% Sodium Azide