

Product Datasheet

Campylobacter Antibody HM806

Product Name	Campylobacter Antibody
Catalog Number	HM806
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% by SDS-PAGE
Buffer	Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with <i>H. pylori</i> , <i>E. coli</i> , <i>S. enteritidis</i> , <i>S. paratyphi</i> , <i>S. typhi</i> , <i>S. typhimurium</i> , <i>S. boydii</i> , <i>S. dysenteriae</i> , <i>S. flexneri</i> , <i>S. sonnei</i> , <i>Y. enterocolitica</i> , <i>L. monocytogenes</i> , <i>S. aureus</i>
Matched Pairs	HM805 (detection)
Preservatives	0.09% Sodium Azide