

Product Datasheet

Cryptosporidium Antibody HM807

Product Name	Cryptosporidium Antibody
Catalog Number	HM807
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified Native <i>C. parvum</i>
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with <i>E. dispar</i> , <i>E. hystolitica</i> , <i>G. lamblia</i>
Specificity	Cryptosporidium
Matched Pairs	HM807 (self-pairs)
Preservatives	0.09% Sodium Azide