

Product Datasheet

Human IL-17A/F ELISA Kit 850.990

Product Name	Human IL-17A/F ELISA Kit
Catalog Number	850.990
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	15.6-500 pg/ml
Sensitivity	6.6 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reaction with IL-17A, IL-17B, IL-17D, IL-17E, IFN-g, Gal-1 and IL-33. Shows cross reactivity with IL-17F.
Specificity	Recognizes both natural and recombinant human IL-17A/F
Sample Types	Serum, Plasma, Cell culture supernatant