

Product Datasheet

Entamoeba Antibody HM808

Product Name	Entamoeba Antibody
Catalog Number	HM808
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with <i>C. parvum</i> , <i>G. lamblia</i>
Matched Pairs	HM809 (capture)
Preservatives	0.09% Sodium Azide