

Product Datasheet

Lactoferrin Antibody HM811

Product Name	Lactoferrin Antibody
Catalog Number	HM811
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with bovine lactoferrin, human/pig/bovine hemoglobin, human calprotectin, human/pig/bovine transferrin
Matched Pairs	HM810 (detection)
Preservatives	0.09% Sodium Azide