

## Product Datasheet

### Lactoferrin Antibody HM811

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Lactoferrin Antibody   |
| <b>Catalog Number</b>      | HM811  |
| <b>Species/Host</b>        | Mouse  |
| <b>Tested Applications</b> | LF, ELISA  |
| <b>Brand</b>               | EastCoast Bio  |
| <b>Form/Appearance</b>     | Purified Monoclonal Antibody   |
| <b>Storage</b>             | 2-8°C  |
| <b>Isotype</b>             | IgG1   |
| <b>Clonality</b>           | Monoclonal   |
| <b>Purity</b>              | >95% by SDS-PAGE   |
| <b>Buffer</b>              | 10 mM Phosphate Buffered Saline, pH 7.4  |
| <b>Cross Reactivity</b>    | No reactivity with bovine lactoferrin, human/pig/bovine hemoglobin, human calprotectin, human/pig/bovine transferrin |
| <b>Matched Pairs</b>       | HM810 (detection)  |
| <b>Preservatives</b>       | 0.09% Sodium Azide   |