

Product Datasheet

Hemoglobin A1c (HbA1c) Antibody HM813

Product Name	Hemoglobin A1c (HbA1c) Antibody
Catalog Number	HM813
Species/Host	Mouse
Tested Applications	LF, ELISA, WB, LTIA
Immunogen	Native HbA1c with peptides
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM814 (capture)
Preservatives	0.05% Proclin 300