

Product Datasheet

Brain Natriuretic Peptide (BNP) Antibody HM816

Product Name	Brain Natriuretic Peptide (BNP) Antibody
Catalog Number	HM816
Species/Host	Mouse
Tested Applications	LF
Immunogen	Recombinant BNP
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM817 (capture)
Preservatives	0.05% Proclin 300