

Product Datasheet

Brain Natriuretic Peptide (BNP) Antibody HM818

Product Name	Brain Natriuretic Peptide (BNP) Antibody
Catalog Number	HM818
Species/Host	Mouse
Tested Applications	ELISA, CLIA
Immunogen	Recombinant BNP
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Epitope	26-32 a.a.
Purity	>90%
Buffer	Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM819 (capture)
Preservatives	0.05% Proclin 300