

## Product Datasheet

### Brain Natriuretic Peptide (BNP) Antibody HM819

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Brain Natriuretic Peptide (BNP) Antibody                                  |
| <b>Catalog Number</b>      | HM819   |
| <b>Species/Host</b>        | Mouse   |
| <b>Tested Applications</b> | ELISA, CLIA   |
| <b>Immunogen</b>           | Recombinant BNP   |
| <b>Brand</b>               | EastCoast Bio   |
| <b>Form/Appearance</b>     | Purified Monoclonal Antibody  |
| <b>Storage</b>             | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| <b>Isotype</b>             | IgG2a   |
| <b>Clonality</b>           | Monoclonal  |
| <b>Epitope</b>             | 5-13 a.a.   |
| <b>Purity</b>              | >90%  |
| <b>Buffer</b>              | Phosphate Buffered Saline, pH 7.4   |
| <b>Matched Pairs</b>       | HM818 (detection)   |
| <b>Preservatives</b>       | 0.05% Proclin 300   |