

Product Datasheet

Lipoprotein-associated Phospholipase A2 (Lp-PLA2) Antibody HM821

Product Name	Lipoprotein-associated Phospholipase A2 (Lp-PLA2) Antibody
Catalog Number	HM821
Species/Host	Mouse
Tested Applications	LF, ELISA, WB
Immunogen	Native Lp-PLA2
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Lp-PLA2
Matched Pairs	HM820 (detection)
Preservatives	0.05% Proclin 300