

Product Datasheet

Vitamin B12 Antibody HM831

Product Name	Vitamin B12 Antibody
Catalog Number	HM831
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Vitamin B12-BSA
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	LA484 (capture)
Preservatives	0.05% Sodium Azide