

## Product Datasheet

### Serum Amyloid A (SAA), Horse Antibody HM840

<b>Product Name</b>	Serum Amyloid A (SAA), Horse Antibody
<b>Catalog Number</b>	HM840
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA, WB
<b>Immunogen</b>	Horse SAA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term (<2 weeks): 2-8°C. Long Term: -20°C
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>90%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Matched Pairs</b>	HM840 (self-pairs)
<b>Preservatives</b>	0.1% Sodium Azide