

Product Datasheet

Serum Amyloid A (SAA), Horse Antibody HM840

Product Name	Serum Amyloid A (SAA), Horse Antibody
Catalog Number	HM840
Species/Host	Mouse
Tested Applications	LF, ELISA, WB
Immunogen	Horse SAA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (<2 weeks): 2-8°C. Long Term: -20°C
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM840 (self-pairs)
Preservatives	0.1% Sodium Azide