

## Product Datasheet

### Dengue Virus (DENV-pan NS1) Antibody HM843

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Dengue Virus (DENV-pan NS1) Antibody                                      |
| <b>Catalog Number</b>      | HM843   |
| <b>Species/Host</b>        | Mouse   |
| <b>Tested Applications</b> | LF, ELISA   |
| <b>Immunogen</b>           | Purified Dengue Virus   |
| <b>Brand</b>               | EastCoast Bio   |
| <b>Form/Appearance</b>     | Purified Monoclonal Antibody  |
| <b>Storage</b>             | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| <b>Isotype</b>             | IgG1  |
| <b>Clonality</b>           | Monoclonal  |
| <b>Purity</b>              | >90%  |
| <b>Buffer</b>              | 0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride             |
| <b>Specificity</b>         | NS1 protein of all 4 serotypes  |
| <b>Matched Pairs</b>       | HM349 or HM388 (capture)  |
| <b>Preservatives</b>       | 0.05% Sodium Azide  |