

## Product Datasheet

### Bovine Growth Hormone (rBGH) Antibody HM844

<b>Product Name</b>	Bovine Growth Hormone (rBGH) Antibody
<b>Catalog Number</b>	HM844
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA
<b>Immunogen</b>	rBGH
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Ion Exchange Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>90% by SDS-PAGE
<b>Buffer</b>	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
<b>Cross Reactivity</b>	Pituitary Bovine Growth Hormone (pBGH), Native Bovine Growth Hormone (nBGH)
<b>Specificity</b>	Recombinant Bovine Growth Hormone
<b>Matched Pairs</b>	HM845 or HM847 (detection)
<b>Preservatives</b>	0.05% Sodium Azide