

Product Datasheet

Bovine Growth Hormone (rBGH) Antibody HM845

Product Name	Bovine Growth Hormone (rBGH) Antibody
Catalog Number	HM845
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	rBGH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	Pituitary Bovine Growth Hormone (pBGH), Native Bovine Growth Hormone (nBGH)
Specificity	Recombinant Bovine Growth Hormone
Matched Pairs	HM844 or HM846 (capture)
Preservatives	0.05% Sodium Azide