

Product Datasheet

Bovine Growth Hormone (rBGH) Antibody HM847

| | |
|----------------------------|---|
| Product Name | Bovine Growth Hormone (rBGH) Antibody |
| Catalog Number | HM847 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | rBGH |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | >90% by SDS-PAGE |
| Buffer | 0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride |
| Cross Reactivity | Pituitary Bovine Growth Hormone (pBGH), Native Bovine Growth Hormone (nBGH) |
| Specificity | Recombinant Bovine Growth Hormone |
| Matched Pairs | HM844 or HM846 (capture) |
| Preservatives | 0.05% Sodium Azide |