

Product Datasheet

Creatine Kinase MB (CK-MB) Antibody HM848

Product Name	Creatine Kinase MB (CK-MB) Antibody
Catalog Number	HM848
Species/Host	Mouse
Tested Applications	LF, ELISA, LTIA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM849 (detection)
Preservatives	0.1% Proclin 300