

Product Datasheet

Parvovirus, Feline Antibody HM853

Product Name	Parvovirus, Feline Antibody
Catalog Number	HM853
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Native Feline Parvovirus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Matched Pairs	HM852 (detection)
Preservatives	0.05% Sodium Azide