

## Product Datasheet

### Parvovirus, Canine Antibody HM860

<b>Product Name</b>	Parvovirus, Canine Antibody
<b>Catalog Number</b>	HM860
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA, WB
<b>Immunogen</b>	Canine Parvovirus
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term ( $\geq 2$ weeks) 2-8°C. Long Term: -20°C or below
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Matched Pairs</b>	HM861 (detection)
<b>Preservatives</b>	0.09% Sodium Azide