

Product Datasheet

Parvovirus, Canine Antibody HM861

Product Name	Parvovirus, Canine Antibody
Catalog Number	HM861
Species/Host	Mouse
Tested Applications	LF, ELISA, WB
Immunogen	Canine Parvovirus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≥ 2 weeks) 2-8°C. Long Term: -20°C or below
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM860 (capture)
Preservatives	0.09% Sodium Azide