

Product Datasheet

Serum Amyloid A (SAA), Feline Antibody HM866

Product Name	Serum Amyloid A (SAA), Feline Antibody
Catalog Number	HM866
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM864 or HM865 (capture)