

## Product Datasheet

### Pepsinogen II (PGII) Antibody HM876

<b>Product Name</b>	Pepsinogen II (PGII) Antibody
<b>Catalog Number</b>	HM876
<b>Species/Host</b>	CHO Cells
<b>Tested Applications</b>	ELISA, CLIA, FIA, TIA
<b>Immunogen</b>	Recombinant pepsinogen II
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>90%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Specificity</b>	PGII
<b>Matched Pairs</b>	HM877 (detection)
<b>Preservatives</b>	0.1% Proclin 300