

Product Datasheet

Pepsinogen II (PGII) Antibody HM877

Product Name	Pepsinogen II (PGII) Antibody
Catalog Number	HM877
Species/Host	CHO Cells
Tested Applications	ELISA, CLIA, FIA, TIA
Immunogen	Recombinant pepsinogen II
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 6.0
Specificity	PGII
Matched Pairs	HM876 (capture)
Preservatives	0.1% Proclin 300