

Product Datasheet

Rift Valley Fever (RVFV-NP) Antibody HM878

Product Name	Rift Valley Fever (RVFV-NP) Antibody
Catalog Number	HM878
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant RVFV
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Specificity	RVFV nucleoprotein
Matched Pairs	HM879 (detection)
Preservatives	0.05% Sodium Azide