

Product Datasheet

Pepsinogen I (PGI) Antibody HM887

Product Name	Pepsinogen I (PGI) Antibody
Catalog Number	HM887
Species/Host	CHO Cells
Tested Applications	ELISA, CLIA, FIA, TIA
Immunogen	Recombinant pepsinogen I
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffer with 150 mM NaCl, pH 6.0
Specificity	PGI
Matched Pairs	HM886 (detection)
Preservatives	0.1% Proclin 300