

Product Datasheet

Anti-h Trypsinogen-2 8607 SPTN-5 100108

Product Name	Anti-h Trypsinogen-2 8607 SPTN-5
Catalog Number	100108
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	1.6×10^8 1/M
Buffer	50 mM Na-citrate, pH 7.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.0-6.9
Cross Reactivity	Human Trypsin-1 2 %, Human Trypsin-2 6 %
Specificity	Antibody recognizes human trypsinogen-2
Shelf Life	12 months