

Product Datasheet

Troponin I (85-92) Antibody HM898

Product Name	Troponin I (85-92) Antibody
Catalog Number	HM898
Species/Host	CHO Cells
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Epitope	85-92 a.a.
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Troponin I
Matched Pairs	HM897 (capture) and HM896 (detection)
Preservatives	0.1% Proclin 300