

Product Datasheet

NT-proBNP (43-50) Antibody HM899

Product Name	NT-proBNP (43-50) Antibody
Catalog Number	HM899
Species/Host	Sheep
Tested Applications	LF, CLIA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10mM Phosphate Buffered Saline, pH 7.4
Specificity	NT-proBNP (a.a. 43-50)
Matched Pairs	HM148 or HM210 (detection)
Preservatives	0.1% Proclin 300