

Product Datasheet

Senecavirus A (SVA) Antibody HM903

Product Name	Senecavirus A (SVA) Antibody
Catalog Number	HM903
Species/Host	Rat
Tested Applications	LF, ELISA
Immunogen	Recombinant Senecavirus A
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.15M Potassium Phosphate, 0.85% NaCl, pH 7.2
Matched Pairs	HM900 (capture)
Preservatives	0.05% Sodium Azide