

Product Datasheet

CYFRA21-1 (Cytokeratin Fragment) Antibody HM907

Product Name	CYFRA21-1 (Cytokeratin Fragment) Antibody
Catalog Number	HM907
Species/Host	Mouse
Tested Applications	ELISA, CLIA
Immunogen	Recombinant CYFRA21-1
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffer, pH 7.4 with 150 mM NaCl
Matched Pairs	HM908 (detection)
Preservatives	0.1% Proclin 300