

Product Datasheet

Inhibin B (INHB) Antibody HM909

Product Name	Inhibin B (INHB) Antibody
Catalog Number	HM909
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Recombinant INHB
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: ≤-20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffer, pH 7.4 with 150 mM NaCl
Specificity	Inhibin B β subunit
Matched Pairs	HM910 (detection)
Preservatives	0.1% Proclin 300