

Product Datasheet

Inhibin B (INHB) Antibody HM910

Product Name	Inhibin B (INHB) Antibody
Catalog Number	HM910
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Recombinant INHB
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	-20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffer, pH 7.4 with 150 mM NaCl
Specificity	Inhibin B α subunit
Matched Pairs	HM909 (capture)
Preservatives	0.1% Proclin 300