

Product Datasheet

Biotin Antibody HM911

Product Name	Biotin Antibody
Catalog Number	HM911
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Biotin-KLH
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	Phosphate Buffered Saline
Specificity	Biotin
Preservatives	0.09% Sodium Azide