

## Product Datasheet

### Human Immunoglobulin E (IgE) Antibody HM912

<b>Product Name</b>	Human Immunoglobulin E (IgE) Antibody
<b>Catalog Number</b>	HM912
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Purified Human IgE
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Cross Reactivity</b>	No reactivity with Human IgA, IgD, IgG, IgM
<b>Specificity</b>	Human IgE
<b>Preservatives</b>	0.09% Sodium Azide