

Product Datasheet

Human Immunoglobulin E (IgE) Antibody HM912

Product Name	Human Immunoglobulin E (IgE) Antibody
Catalog Number	HM912
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Purified Human IgE
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Human IgA, IgD, IgG, IgM
Specificity	Human IgE
Preservatives	0.09% Sodium Azide